

ELECTRO
PORATION



Complete with the very best of professional cosmetics

An entire line of cosmetics designed by Wavemed offers you everything you need to meet your client's needs. To restore tone and freshness to the skin of the face and body, Wavemed offers four products specifically designed to be deeply transmitted by the electroporator. The device enhances the effectiveness of the products.

ACIDO IALURONICO

face and body

Wavemed Hyaluronic acid

is an excellent remedy against wrinkles and skin dehydration. By stimulating cellular capabilities, it reactivates mature skin, restoring elasticity and softness.

VITAMINA B3

face and body

Vitamin B3

counteracts the formation of skin spots and, thanks to its antioxidant properties, it has a toning and smoothing effect.

EDERA E CENTELLA

face and body

Ivy and Centella

are plant extracts which have a firming and draining effect on tissues. Besides giving tone to the skin the extracts counteract venous insufficiency and leg swelling.

FOSFATIDILCOLINA

face and body

Phosphatidylcholine

reduces localized fat favouring the emulsification of fats and their natural expulsion. It is widely used in anti-cellulite treatments, but it is also a great way to restore vigour and tone to the skin.



The hand-piece that guarantees results

To best counter the blemishes of ageing, the active ingredients of the cosmetic products are applied by passing them through the skin. To ensure the effectiveness of the treatment, the hand-piece does not have any metal components. This ensures that the deep absorption of the active ingredients is due to the polarization of the molecules of the cosmetic product. The results are the final confirmation.

The optimal combination

Thanks to the modular design of Modula, the electroporator may be paired with the radio frequency module, obtaining the ideal platform for the rejuvenation of face and body skin. The heat transmitted in depth by the Rf favours the passage of molecules into the dermis, accelerating the benefits. You are able to benefit from the two technologies with the same hand-piece. It is the only electroporator with double effectiveness!



Modula

ELECTROPORATION

General

Display Touch Screen 10,4''

Power source	110 – 240 VAC 50-60 Hz
Absorption	150 VA
Device type	BF
Medical Device	Class IIb - CE 1936
Insulation class	I
Fuses	2 x 250V T4A
Casing	Fire proof ABS
Frequency	800-1200 Hz
Current	20 mA
Hand-piece	Roll-on Model
Connectivity	Wifi, LAN, USB
Size l x h x d	400 x 310 x 230 mm
Accessories	Cart, trolley
Weight	8 kg
Accuracy of displayed values	± 20 %